




ENDEAVOUR BARIATRIC TOILET SURROUND



SKU: E208

PRODUCT DESCRIPTION

Purpose & Benefits

The **Endeavour Bariatric Toilet Surround** provides sturdy and stable support for users requiring assistance when sitting or standing from the toilet. Its strong, adjustable frame offers added safety and confidence in daily bathroom use.

Features

- Provides firm support when transferring on and off the toilet
- Height adjustable to suit user preference and toilet height
- Strong and stable frame suitable for bariatric users
- Designed for use in domestic or healthcare environments

Specifications

- **Width:** 600 mm
- **Height Adjustment Range:** 680 – 810 mm
- **Safe Working Load (SWL):** 175 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au