


EQUAGEL GENERAL CUSHION



SKU: CE2-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The EquaGel General is an everyday pressure reducing cushion that is made of a dry-polymer gel that feels almost rubbery to the touch. This gel cushion is constructed in a grid system and when a user sits on the cushion creates a supportive and comfortable surface, aiding pressure relief. The air channels promote air flow and assists in distribution of pressure over a larger surface area. This, along with the breathable cover, optimises skin micro-climate. All EquaGel cushions come with breathable, water-resistant and multi-stretch covers. Featuring Velcro straps for wheelchair attachment, these cushions are long-lasting, low maintenance and easy to use. EquaGel cushions are maintenance-free, do not require set up and are flexible for use on slung wheelchairs seats.

KEY FEATURES

- Ideal for mild to moderate cushioning needs
- Contoured pelvic area for equalising pressure redistribution
- Waterproof, flame retardant, vapour-permeable cover with hook-and-loop sections on base for wheelchair security
- 50 mm high, with 32 mm of buckling gel under pressure points
- No valve or pump. Zero maintenance
- Vapour permeable, waterproof cover. Maintains dry skin.
- Offer clinical advantages over inflatable and foam cushions.
- 12+ degrees cooler than inflatable cushions.
- Foam cushions amplify vibrations, EquaGel absorbs. This has a big impact on skin and tissue

integrity for those using scooters.

- All components are machine washable.

SPECIFICATIONS

- 1.75" total height
- Safe Working Load - Unlimited

Buy this EquaGel cushion online at Medimart Australia using [NDIS payment method](#), AfterPay or PayPal.

[View Product Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

