



## **EQUAGEL PROTECTOR CUSHION**

**\$450.00** Sale Available Online Only!



**SKU:** CE21-H

## **PRODUCT DESCRIPTION**

### **EQUAGEL PROTECTOR CUSHION - PURPOSE AND BENEFITS**

The EquaGel Protector Cushion is a cutting-edge seating solution designed to provide optimum comfort and pressure redistribution for individuals who spend extended periods sitting. Combining gel and column-buckling technologies, this cushion is particularly favoured among wheelchair users and others in need of advanced pressure relief.

### **FEATURES**

- **Gel-Column Technology:** The cushion uses a gel matrix of columns that buckle and distribute weight evenly, effectively reducing peak pressures.
- **Deep Immersion:** Offers superior deep tissue protection, which aids in reducing pressure on the seated area.
- **Two-Layer Design:** Features a thin, soft top layer and a thicker, firmer bottom layer for maximum comfort and support.
- **Fluid-Resistant Cover:** Comes with a waterproof and easy-to-clean cover, enhancing the cushion's durability.
- **Maintenance-Free:** Requires no pumps, adjustments, or regular maintenance.
- **Stability Trench:** The design aids in maintaining pelvic positioning, ensuring user stability and postural support.

### **BENEFITS OF THE EQUAGEL PROTECTOR CUSHION**

- **Pressure Relief:** Helps in preventing and treating pressure sores and ulcers.

- **Comfort:** Provides prolonged seating comfort, especially for those with limited mobility.
- **Durability:** Made of high-quality materials, ensuring a long lifespan even with daily use.
- **Therapeutic:** Ideal for individuals recovering from surgeries or facing conditions that affect their sitting posture.

## **APPLICATIONS**

The EquaGel Protector Cushion is designed for individuals who:

- Use wheelchairs regularly.
- Have a history of or are at risk of developing pressure sores.
- Need added seating comfort due to various health conditions.
- Spend extended hours seated, either due to professional requirements or mobility challenges.

## **ADDITIONAL INFORMATION**

The EquaGel Protector Cushion is 2.5 inches thick and has a contoured pelvic area to more evenly distribute weight. It has a fluid-resistant, vapor-permeable, and flame-resistant cover. The cushion also has a velcro bottom for wheelchair users. The cushion also has a safe working load of 270 kilograms, making it safe for all users needing the support.