

EQUAGEL PROTECTOR CUSHION



SKU: CE21-H



PRODUCT DESCRIPTION

EQUAGEL PROTECTOR CUSHION - PURPOSE AND BENEFITS

The EquaGel Protector Cushion is a cutting-edge seating solution designed to provide optimum comfort and pressure redistribution for individuals who spend extended periods sitting. Combining gel and column-buckling technologies, this cushion is particularly favoured among wheelchair users and others in need of advanced pressure relief.

FEATURES

- **Gel-Column Technology:** The cushion uses a gel matrix of columns that buckle and distribute weight evenly, effectively reducing peak pressures.
- **Deep Immersion:** Offers superior deep tissue protection, which aids in reducing pressure on the seated area.
- **Two-Layer Design:** Features a thin, soft top layer and a thicker, firmer bottom layer for maximum comfort and support.
- **Fluid-Resistant Cover:** Comes with a waterproof and easy-to-clean cover, enhancing the cushion's durability.
- **Maintenance-Free:** Requires no pumps, adjustments, or regular maintenance.
- **Stability Trench:** The design aids in maintaining pelvic positioning, ensuring user stability and postural support.

BENEFITS OF THE EQUAGEL PROTECTOR CUSHION

- **Pressure Relief:** Helps in preventing and treating pressure sores and ulcers.

- **Comfort:** Provides prolonged seating comfort, especially for those with limited mobility.
- **Durability:** Made of high-quality materials, ensuring a long lifespan even with daily use.
- **Therapeutic:** Ideal for individuals recovering from surgeries or facing conditions that affect their sitting posture.

APPLICATIONS

The EquaGel Protector Cushion is designed for individuals who:

- Use wheelchairs regularly.
- Have a history of or are at risk of developing pressure sores.
- Need added seating comfort due to various health conditions.
- Spend extended hours seated, either due to professional requirements or mobility challenges.

ADDITIONAL INFORMATION

The EquaGel Protector Cushion is 2.5 inches thick and has a contoured pelvic area to more evenly distribute weight. It has a fluid-resistant, vapor-permeable, and flame-resistant cover. The cushion also has a velcro bottom for wheelchair users. The cushion also has a safe working load of 270 kilograms, making it safe for all users needing the support.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au