

## ETAC AKTIV REACHER WITH HOOK

**\$80.00**



**SKU:** 80505006

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

The Etac Aktiv reacher with hook is a convenient assistive tool that proves invaluable for dressing and undressing tasks, thanks to its pull and push capabilities.

## FEATURES

- This reacher shares the same lightweight construction and excellent functionality as other reachers on the market.
- The handles are thoughtfully designed to leverage the strength of the hand, offering an effective grip for users.
- The jaws of the reacher are coated with durable rubber, ensuring ample friction to minimize the required amount of strength while providing a secure hold on objects.

## SPECIFICATIONS

Item Length	Item number	Bracket material	Design	Frame material	Handle material	Jaw material	Lifting capacity, max (kg)	Weight (g)
45 cm	80505006	Polyamide	Ergonomidesign	Aluminium	Jaw bracket: Polyamide	Polycarbonate/ABS	2	117