

## ETAC AKTIV REACHER WITH POWER GRIP AND HOOK



**SKU:** 80505005

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Aktiv reacher with Power Grip and Hook is specifically designed to provide support for individuals with weak hands and limited lifting ability. It features a forearm support that assists in maintaining stability while using the reacher.

### FEATURES

- Similar to other Aktiv reachers, the jaws of this model are turnable, allowing for versatile use.
- Equipped with hooks that enable both pushing and pulling of items
- To enhance durability and improve gripping capabilities, the jaws are coated with hardwearing rubber.
- A magnet is incorporated into the reacher's design, enabling it to pick up small metallic items conveniently.

### SPECIFICATIONS

Item Length	Item number	Design	Frame material	Lifting capacity, max (kg)	Material	Weight (g)
60 cm	80505005	Ergonomidesign	Aluminium	2	Handle and jaws: Polycarbonate/ABS. Forearm support and jaw bracket: Polyamide	139

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)