

ETAC AKTIV REACHER WITH POWER GRIP AND HOOK

\$73.00



SKU: 80505005



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Aktiv reacher with Power Grip and Hook is specifically designed to provide support for individuals with weak hands and limited lifting ability. It features a forearm support that assists in maintaining stability while using the reacher.

FEATURES

- Similar to other Aktiv reachers, the jaws of this model are turnable, allowing for versatile use.
- Equipped with hooks that enable both pushing and pulling of items
- To enhance durability and improve gripping capabilities, the jaws are coated with hardwearing rubber.
- A magnet is incorporated into the reacher's design, enabling it to pick up small metallic items conveniently.

SPECIFICATIONS

Item Length	Item number	Design	Frame material	Lifting capacity, max (kg)	Material	Weight (g)
60 cm	80505005	Ergonomidesign	Aluminium	2	Handle and jaws: Polycarbonate/ABS. Forearm support and jaw bracket: Polyamide	139