

## ETAC AKTIV REACHER WITH POWER GRIP AND HOOK

**\$73.00**



**SKU:** 80505005

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Aktiv reacher with Power Grip and Hook is specifically designed to provide support for individuals with weak hands and limited lifting ability. It features a forearm support that assists in maintaining stability while using the reacher.

### FEATURES

- Similar to other Aktiv reachers, the jaws of this model are turnable, allowing for versatile use.
- Equipped with hooks that enable both pushing and pulling of items
- To enhance durability and improve gripping capabilities, the jaws are coated with hardwearing rubber.
- A magnet is incorporated into the reacher's design, enabling it to pick up small metallic items conveniently.

### SPECIFICATIONS

Item Length	Item number	Design	Frame material	Lifting capacity, max (kg)	Material	Weight (g)
60 cm	80505005	Ergonomidesign	Aluminium	2	Handle and jaws: Polycarbonate/ABS. Forearm support and jaw bracket: Polyamide	139