

ETAC AKTIV REACHER

\$60.00 - \$100.00 Price range: \$60.00 through \$100.00



SKU: 8050500-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The handles of the Etac Aktiv reacher are thoughtfully designed to optimize the utilization of hand strength and provide an effective grip.

These handles are slightly angled, allowing for easier maneuverability, and require only a small hand grip to operate. Additionally, they feature a support for the index finger, which enhances precision and control during use.

FEATURES

- Handles are slightly angled, allowing for easier maneuverability
- Features a support for the index finger, which enhances precision and control during use.
- To enhance durability and improve gripping capabilities, the jaws are coated with hardwearing rubber.

SPECIFICATIONS

Item Length	Item number	Bracket material	Design	Frame material	Handle material	Jaw material	Lifting capacity, max (kg)	Weight (g)
45 cm	80505004	Polyamide	Ergonomidesign	Aluminium	Polycarbonate/ABS	Polycarbonate/ABS	2	111
60 cm	80505007	Polyamide	Ergonomidesign	Aluminium	Polycarbonate/ABS	Polycarbonate/ABS	2	120
70 cm	80505003	Polyamide	Ergonomidesign	Aluminium	Polycarbonate/ABS	Polycarbonate/ABS	2	128
80 cm	80505008	Polyamide	Ergonomidesign	Aluminium	Polycarbonate/ABS	Polycarbonate/ABS	2	138
90 cm	80505009	Polyamide	Ergonomidesign	Aluminium	Polycarbonate/ABS	Polycarbonate/ABS	2	150

