



## ETAC CLEAN MOBILE SHOWER COMMODE - HEIGHT ADJUSTABLE



**SKU:** 80229208

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

Etac Clean Height Adjustable shower commode chair offers the convenience of adjustable seat height, allowing for a customized fit to meet the user's specific needs. This adjustment capability makes transfers easier and ensures a perfect fit over the toilet. The height adjustment process is simple and tool-free.

### FEATURES

- The Etac Clean shower commode chair is specifically designed to fit perfectly over wall-mounted toilets, thanks to its upright operating handle. This feature ensures seamless integration with various toilet configurations, allowing for easy and convenient use of the chair during toileting tasks.
- Etac Clean has a smart design that provides a roomy and comfortable sitting area while maintaining a slim profile.
- The foot support of the Etac Clean shower commode chair is designed with a soft, rounded shape, providing both comfort and protection to the user's feet.
- The Etac Clean shower commode chair has an open back for easy access during hygiene activities.

### SPECIFICATIONS

Item Name	Item number	Colour	Maximum user weight (kg)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)
-----------	-------------	--------	--------------------------	-----------------	------------------	-----------------	-------------	------------

Clean height-adjustable, 4 lockable wheels	80229208	White	130	440	475 - 600	480	14.2	540
--	----------	-------	-----	-----	-----------	-----	------	-----

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)