



ETAC CLEAN MOBILE SHOWER COMMODE



SKU: 80229-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Clean shower commode chair is designed for easy operation, made possible by the width of the push handle. The chair's clever design is evident throughout its entirety. The upright handle ensures a perfect fit, even over wall-mounted toilets, while the chair's open back enhances access and facilitates use.

The Clean shower commode chair features a unique foot support that offers both comfort and practicality. The foot support can be conveniently pushed under the seat when not in use. The soft curved surface of the foot support ensures comfortable foot placement, and the side protection of the foot support keeps the feet well-protected.

Additional Functionality in Other Variants All other variants of the Clean shower commode chair retain these features while also incorporating additional functionality to cater to specific needs.

FEATURES

- Etac Clean has a smart design that provides a roomy and comfortable sitting area while maintaining a slim profile.
- The foot support of the Etac Clean shower commode chair is designed with a soft, rounded shape, providing both comfort and protection to the user's feet.
- The Etac Clean shower commode chair is designed to fit effortlessly over wall-mounted toilets, thanks to its upright operating handle.
- The Etac Clean shower commode chair features an open seat design for easy access and hygiene purposes, available with front or back openings.

SPECIFICATIONS

Item Name	Item number	Colour	Maximum user weight (kg)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Width (mm)
Clean 49 cm, 4 lockable castors	80229298	White	130	440	490	480	500
Clean 55 cm, 4 lockable castors	80229271	White	130	440	550	480	500
Clean 49 cm, 4 lockable castors	80229210	Lagoon green	130	440	490	480	500
Clean 55 cm, 4 lockable castors	80229211	Lagoon green	130	440	550	480	500
Clean 55 cm, front-opened seat, 4 lockable castors	80229321	Lagoon green	130	440	550	480	500
Clean 49 cm, front-opened seat, 4 lockable castors	80229325	Lagoon green	130	440	490	480	500

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au