




## **ETAC CLEAN MOBILE SHOWER COMMODE**

**\$1,100.00 - \$1,600.00**



**SKU: 80229-H**

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Etac Clean shower commode chair is designed for easy operation, made possible by the width of the push handle. The chair's clever design is evident throughout its entirety. The upright handle ensures a perfect fit, even over wall-mounted toilets, while the chair's open back enhances access and facilitates use.

The Clean shower commode chair features a unique foot support that offers both comfort and practicality. The foot support can be conveniently pushed under the seat when not in use. The soft curved surface of the foot support ensures comfortable foot placement, and the side protection of the foot support keeps the feet well-protected.

**Additional Functionality in Other Variants** All other variants of the Clean shower commode chair retain these features while also incorporating additional functionality to cater to specific needs.

### **FEATURES**

- Etac Clean has a smart design that provides a roomy and comfortable sitting area while maintaining a slim profile.
- The foot support of the Etac Clean shower commode chair is designed with a soft, rounded shape, providing both comfort and protection to the user's feet.
- The Etac Clean shower commode chair is designed to fit effortlessly over wall-mounted toilets, thanks to its upright operating handle.
- The Etac Clean shower commode chair features an open seat design for easy access and hygiene purposes, available with front or back openings.

## SPECIFICATIONS

Item Name	Item number	Colour	Maximum user weight (kg)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Width (mm)
Clean 49 cm, 4 lockable castors	80229298	White	130	440	490	480	500
Clean 55 cm, 4 lockable castors	80229271	White	130	440	550	480	500
Clean 49 cm, 4 lockable castors	80229210	Lagoon green	130	440	490	480	500
Clean 55 cm, 4 lockable castors	80229211	Lagoon green	130	440	550	480	500
Clean 55 cm, front-opened seat, 4 lockable castors	80229321	Lagoon green	130	440	550	480	500
Clean 49 cm, front-opened seat, 4 lockable castors	80229325	Lagoon green	130	440	490	480	500