




## ETAC CLEAN SELF PROPELLED MOBILE SHOWER COMMODE

**\$2,000.00 - \$2,306.00** Price range: \$2,000.00 through \$2,306.00



**SKU:** 802292-H

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

The Etac Clean shower commode chair is designed with 24-inch wheels strategically positioned for effortless maneuverability while ensuring user safety and comfort. Transferring to and from the chair is made easier by the presence of a large, unobstructed seat in front of the wheels. The placement of the brake handles below the seat allows for unimpeded transfers and remains easily accessible. This thoughtful design consideration ensures smooth and convenient usage of the chair while maintaining user comfort and safety.

## FEATURES

- Self-propelled making it easy to operate for the independent user.
- The open seat of the Etac Clean shower commode chair allows for easy access during hygiene routines, with the option of front or back openings.
- The foot support of the Etac Clean shower commode chair has a soft, rounded shape for extra comfort and protection.
- The flat seat and brake position of the Etac Clean shower commode chair facilitate safe and simple lateral transfers.

## SPECIFICATIONS

Item Name	Item number	Colour	Maximum user weight (kg)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)

Clean 55 cm with 24" rear wheels	80229276	White	130	440	550	480	20	670
Clean 55 cm with 24" rear wheels	80229216	Lagoon green	130	440	550	480	20	670