



## ETAC CLEAN SELF PROPELLED MOBILE SHOWER COMMODE



**SKU:** 802292-H

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Clean shower commode chair is designed with 24-inch wheels strategically positioned for effortless maneuverability while ensuring user safety and comfort. Transferring to and from the chair is made easier by the presence of a large, unobstructed seat in front of the wheels. The placement of the brake handles below the seat allows for unimpeded transfers and remains easily accessible. This thoughtful design consideration ensures smooth and convenient usage of the chair while maintaining user comfort and safety.

### FEATURES

- Self-propelled making it easy to operate for the independent user.
- The open seat of the Etac Clean shower commode chair allows for easy access during hygiene routines, with the option of front or back openings.
- The foot support of the Etac Clean shower commode chair has a soft, rounded shape for extra comfort and protection.
- The flat seat and brake position of the Etac Clean shower commode chair facilitate safe and simple lateral transfers.

### SPECIFICATIONS

Item Name	Item number	Colour	Maximum user weight (kg)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)

Clean 55 cm with 24" rear wheels	80229276	White	130	440	550	480	20	670
Clean 55 cm with 24" rear wheels	80229216	Lagoon green	130	440	550	480	20	670

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)