



### ETAC CLOO TOILET SEAT RAISER

**\$408.00 - \$412.00** Price range: \$408.00 through \$412.00



**SKU:** 803012-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

Introducing the Etac Cloo, a height-adjustable toilet seat raiser that offers a perfect balance of flexibility and simplicity.

The Cloo can be easily adjusted to three different heights, allowing users to find the most comfortable and suitable position. In addition, it features an option to set a slight anterior tilt, which can be beneficial for individuals with reduced mobility in the hip and knee joints, making sitting down and standing up easier.

Designed with both functionality and aesthetics in mind, the Cloo boasts a discreet and soft design that seamlessly blends into any bathroom interior. This ensures that it not only provides practical assistance but also enhances the overall ambience of the space.

### FEATURES

- The Etac Cloo has a discreet and soft design that effortlessly blends in with the bathroom interior.
- The Etac Cloo allows for easy height adjustment without the need for tools. Spacers can be added or removed effortlessly to achieve the desired seat height.
- The Etac Cloo can be set with a slight anterior tilt to assist individuals with reduced mobility in standing up.
- The Etac Cloo has a removable seat unit with smooth surfaces for easy cleaning. It can be effortlessly detached without the use of tools.

## SPECIFICATIONS

Item Name	Item number	Arm support height from top of seat (mm)	Maximum user weight (kg)	Seat height (mm)	Seat width (mm)	Weight (kg)
Cloo	80301204	-	150	60/100/140*	400	3.3
Cloo with arm supports	80301214	210	150	60/100/140*	400	4.8