



ETAC CLOO TOILET SEAT RAISER



SKU: 803012-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Introducing the Etac Cloo, a height-adjustable toilet seat raiser that offers a perfect balance of flexibility and simplicity.

The Cloo can be easily adjusted to three different heights, allowing users to find the most comfortable and suitable position. In addition, it features an option to set a slight anterior tilt, which can be beneficial for individuals with reduced mobility in the hip and knee joints, making sitting down and standing up easier.

Designed with both functionality and aesthetics in mind, the Cloo boasts a discreet and soft design that seamlessly blends into any bathroom interior. This ensures that it not only provides practical assistance but also enhances the overall ambience of the space.

FEATURES

- The Etac Cloo has a discreet and soft design that effortlessly blends in with the bathroom interior.
- The Etac Cloo allows for easy height adjustment without the need for tools. Spacers can be added or removed effortlessly to achieve the desired seat height.
- The Etac Cloo can be set with a slight anterior tilt to assist individuals with reduced mobility in standing up.
- The Etac Cloo has a removable seat unit with smooth surfaces for easy cleaning. It can be effortlessly detached without the use of tools.

SPECIFICATIONS

Item Name	Item number	Arm support height from top of seat (mm)	Maximum user weight (kg)	Seat height (mm)	Seat width (mm)	Weight (kg)
Cloo	80301204	-	150	60/100/140*	400	3.3
Cloo with arm supports	80301214	210	150	60/100/140*	400	4.8

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au