


## ETAC EASY SHOWER STOOL

**\$200.00 - \$500.00** Price range: \$200.00 through \$500.00



**SKU:** 819010-H



## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Easy Shower Stool is a versatile stool designed to suit a wide range of needs and environments. Its classic rounded seating area ensures easy and comfortable sitting in any position. For added convenience, it is also available with a swivel pad. With its timeless design and adaptability to different interior styles, this stool seamlessly blends into any space, providing both functionality and style. Experience the comfort and convenience of the Etac Easy and transform your living or working environment today.

### FEATURES

- **Comfortable and Safe** - The Etac Easy provides comfortable and safe seating with generously rounded shapes, a non-slip surface, and effective drainage for excess water.
- **Height Adjustable** - The Shower Stool height adjustments are effortless with a click-and-press mechanism and prevents the telescopic legs from turning. Plus, it has no loose parts to worry about misplacing.
- **Easy to assemble** - The assembly of the Etac Easy is tool-free and hassle-free, thanks to the self-locking button that easily clicks into place. The oval shape of the stool prevents the legs from rotating during adjustment, ensuring a secure setup.
- **Stability** - To enhance stability, the legs of the Etac Easy stool are intelligently angled outwards. Additionally, the stool features responsive, non-slip ferrules that adapt to various floor types, ensuring stability even on wet or uneven surfaces.

For a more comfortable seating experience, please check out the [Easy Shower Stool Seat Pad](#)

## SPECIFICATIONS

Item Name	Colour	Material	Maximum user weight (kg)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)
Stool without Swivel	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	420 - 570	400	2.8	495 - 535
Stool with swivel pad	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	420 - 570	400	2.8	495 - 535