


ETAC EDGE SHOWER STOOL



SKU: 818010-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Edge is intelligently designed with a compact, triangular shape, making it a perfect fit for smaller spaces or shower corners. For enhanced stability, users have the option to straddle the seat. Alternatively, they can slightly lower one leg and straddle the seat over that corner, allowing for an open seat angle that can provide relief for aching hips.

FEATURES

- The generous seating areas of the Etac Edge feature softly rounded shapes and a non-slip surface, ensuring both comfort and safety during sitting. Additionally, the stool incorporates special perforations that facilitate the rapid drainage of excess water.
- Adjusting the height of the Etac Edge is a breeze with a simple click-and-press mechanism, facilitated by the oval shape that prevents the telescopic legs from turning during adjustment. Furthermore, there are no loose parts to worry about misplacing.
- The assembly of the Etac Edge requires no tools, as it is designed for effortless setup. The self-locking button securely clicks into place, while the oval shape ensures that the legs remain steady and prevent any turning during adjustment.
- For enhanced stability, the legs of the stools are cleverly angled outwards. Additionally, the stools are equipped with responsive, non-slip ferrules that adapt to different types of flooring, offering exceptional stability even on wet or uneven surfaces.

SPECIFICATIONS

Item Name	Colour	Material	Maximum user weight (kg)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)
Stool	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	420 - 570	450 x 405	2,6	500 - 520
Stool	White	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	420 - 570	450 x 405	2,6	500 - 520
Stool, low	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	345 - 420	450 x 405	2,4	490 - 500

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimart.com.au

