

## ETAC EDGE SHOWER STOOL

**\$185.00 - \$215.00** Price range: \$185.00 through \$215.00



**SKU:** 818010-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Edge is intelligently designed with a compact, triangular shape, making it a perfect fit for smaller spaces or shower corners. For enhanced stability, users have the option to straddle the seat. Alternatively, they can slightly lower one leg and straddle the seat over that corner, allowing for an open seat angle that can provide relief for aching hips.

### FEATURES

- The generous seating areas of the Etac Edge feature softly rounded shapes and a non-slip surface, ensuring both comfort and safety during sitting. Additionally, the stool incorporates special perforations that facilitate the rapid drainage of excess water.
- Adjusting the height of the Etac Edge is a breeze with a simple click-and-press mechanism, facilitated by the oval shape that prevents the telescopic legs from turning during adjustment. Furthermore, there are no loose parts to worry about misplacing.
- The assembly of the Etac Edge requires no tools, as it is designed for effortless setup. The self-locking button securely clicks into place, while the oval shape ensures that the legs remain steady and prevent any turning during adjustment.
- For enhanced stability, the legs of the stools are cleverly angled outwards. Additionally, the stools are equipped with responsive, non-slip ferrules that adapt to different types of flooring, offering exceptional stability even on wet or uneven surfaces.

## SPECIFICATIONS

Item Name	Colour	Material	Maximum user weight (kg)	Seat height (mm)	Seat width (mm)	Weight (kg)	Width (mm)
Stool	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	420 - 570	450 x 405	2,6	500 - 520
Stool	White	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	420 - 570	450 x 405	2,6	500 - 520
Stool, low	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	130	345 - 420	450 x 405	2,4	490 - 500