



ETAC FEED ADJUSTABLE SPOON

\$76.00 Original price was: \$76.00.<u>\$68.00</u>Current price

is: \$68.00.



SKU: 8040330-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The ETAC Feed Adjustable Spoon is specially designed for children to facilitate eating by themselves. Using a power grip, the hand is prevented from sliding thanks to the thumb support. The angle is set with a little screw. The elastic strap is useful for children with poor gripping ability. Etac Feed adjustable spoon is available for both right-handed and left-handed use.

FEATURES

- Lightweight design for ease of use
- Available for left or right-handed users
- Dishwasher safe

SPECIFICATIONS

Weight 40 g **Dimensions** 16 cm

Material White polyamide

Ergonomidesign Design