

ETAC FIX PREPARATION BOARD



SKU: 80501004

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Fix preparation board enables individuals with limited hand strength or single-handed functionality to perform tasks like grating, cutting, whisking, and peeling. The board remains stable with its four friction feet or the suction feet provided. Its jaws securely hold bowls, graters, or food items in place. Additionally, there is a detachable small plate with spikes for holding food pieces. If a flat surface is desired, the plate can be repositioned with the spikes facing downward.

FEATURES

- The jaws of the Fix preparation board offer a reliable grip, ensuring stability while holding a bowl, grater, or food in place during various food preparation tasks.
- The Fix preparation board includes a small plate with spikes specifically designed to securely hold smaller pieces of food. This feature enhances stability and facilitates cutting or chopping tasks. If a flat working surface is desired, the plate can be easily removed and positioned with the spikes pointing downward.
- The clamp of the Fix preparation board provides a secure hold on tins or bowls, offering a convenient and stable platform for food preparation tasks.
- The friction feet of the Fix preparation board effectively prevent any unwanted movement and ensure stability by securely anchoring the board in place during use. This feature is particularly useful on different surfaces, as it minimizes any potential sliding or shifting.
- The Etac Fix Food Preparation Board is crafted using robust and durable materials, making it well-suited for everyday use in the kitchen. Its high-quality construction guarantees long-lasting performance and reliability, ensuring that it can withstand the rigours of regular use over an extended period of time.

SPECIFICATIONS

Item Name	Design	Height (mm)	Length (mm)	Material	Weight (g)	Width (mm)
Etac Fix preparation board	Ergonomidesign. The Fix preparation board is patented	35	315	Chopping board: Polystyrene. Metal parts: Stainless steel. Jaws: PC/ABS	1300	295

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au