


ETAC FLEX GRAB RAIL

\$86.00 - \$212.00Price range: \$86.00 through \$212.00







PURPOSE AND BENEFITS

The Etac Flex grab rail has a sleek, modern design that blends well with contemporary bathrooms. Its modular nature allows for easy customization by connecting multiple rails. The grab rail's shape and material provide a secure grip, even when wet, and the added triangular shape and thumb support enhance safety. The Flex grab rails are available in white and grey.

FEATURES

- You can effortlessly assemble the modular system into a length and shape that perfectly fits your bathroom, offering a tailored solution to your specific needs.
- The grab rail's thoughtfully designed shape and surface ensure a stable grip, even when your hands are wet. This feature enhances safety and confidence during use, minimizing the risk of slipping or losing grip.
- The flexibility of the grab rail system allows you to connect multiple rails according to your
 preferences, enabling you to create a personalized configuration with the desired angles. This
 opens up a range of possibilities, ensuring that the grab rails can be tailored to your specific
 needs and requirements.
- The grab rail offers two attachment options: screws or adhesive. The adhesive option is particularly suitable for individuals who prefer not to or are unable to drill holes in their bathroom walls. This versatile feature allows for easy installation and ensures that the grab rail can be securely mounted without the need for drilling.

SPECIFICATIONS



Item Name	SKU	Colour	Length (mm)	Material	Maximum user weight (kg)
Grab rail, 1x300 mm, screw mount	81707110	White	300	Polypropylene	100
Grab rail, 2x300 mm, screw mount	81707120	White	600	Polypropylene	100
Grab rail, 2x300 mm, screw mount	81707320	Grey	600	Polypropylene	100
Grab rail, 3x300 mm, screw mount	81707330	Grey	900	Polypropylene	100
Extension, screw mount	81706310	Grey	300	Polypropylene	100
Grab rail, 2x300 mm, adhesive mount	81707220	White	600	Polypropylene	100
Extension, adhesive mount	81706210	White	300	Polypropylene	100
Grab rail, 1x300 mm, adhesive mount	81707410	Grey	300	Polypropylene	100
Grab rail, 2x300 mm, adhesive mount	81707420	Grey	600	Polypropylene	100