

ETAC FRESH BATH BOARD



SKU: 816000-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Etac Fresh has established a new benchmark in bath board design with its distinctive features. The wide front edge, positioned on the outer side of the bathtub, offers exceptional support for transfers. In contrast, the mid-section is narrowed off to facilitate easy access for intimate hygiene.

The edge stops are user-friendly and can be easily mounted using the provided plastic coin for secure fastening. Additionally, four friction linings ensure that the bath board remains firmly in place on the bathtub surface. Fresh is available in two different lengths, ensuring a precise fit for various types of bathtubs. It can also be purchased with or without the handle, providing options to suit individual preferences.

FEATURES

- The support handle is designed to provide a comfortable and secure grip, enhancing safety during use. It offers an ergonomic design that promotes ease of handling, allowing users to maintain a stable and confident hold while using the bath board.
- The wide front edge of the bath board, strategically positioned on the outer side of the bathtub, provides exceptional support for transfers. This feature ensures stability and assistance when moving in and out of the bathtub, offering added confidence and security during the transfer process.

SPECIFICATIONS

Item Name	Item number	Adjustability, edge stops (mm)	Design	Length (mm)	Material	Maximum user weight (kg)	Thickness (mm)	Weight (kg)	Width (mm)
-----------	-------------	--------------------------------	--------	-------------	----------	--------------------------	----------------	-------------	------------

Bath board 69 cm	81600011	Max. 636, min. 390	A&E Design	690	Polypropylene. Edge stops: Polypropylene with TPE surface	150	40	2.1	360-270-300
Bath board 69 cm with handle	81600014	Max. 636, min. 390	A&E Design	690	Polypropylene. Edge stops: Polypropylene with TPE surface	150	40	2,1	360-270-300
Bath board 74 cm	81600021	Max. 686, min. 432	A&E Design	740	Polypropylene. Edge stops: Polypropylene with TPE surface	150	40	2.3	360-270-300
Bath board 74 cm with handle	81600024	Max. 686, min. 432	A&E Design	740	Polypropylene. Edge stops: Polypropylene with TPE surface	150	40	2,3	360-270-300

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au