



## ETAC HI-LOO TOILET SEAT RAISER



**SKU:** 8030131-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Hi-Loo with fixed mounting is a discreet toilet seat raiser that is compatible with most toilets. It securely attaches to the toilet, ensuring stability and reliability. The seat itself is ergonomically designed with recesses at the back and front, offering convenience for intimate hygiene purposes.

The arm supports of the Hi-Loo are foldable, but they can also be permanently fixed in the folded-down position if desired. This feature allows for flexibility in limited-space scenarios. Additionally, the Hi-Loo can be used with only one arm support if necessary.

The Hi-Loo fixed model is available in two different heights and can be chosen with or without arm supports, providing options to suit individual preferences and needs.

### FEATURES

- The Etac Hi-Loo Fixed is available in two different heights: 6 cm and 10 cm. This allows users to select the appropriate height that best suits their individual needs and comfort. Whether a slight elevation of 6 cm or a higher raise of 10 cm is desired, the Hi-Loo Fixed provides options for customization and optimal usability.
- The seat unit of the Etac Hi-Loo Fixed is equipped with user-friendly snap locks, making it effortless to detach for easy cleaning.
- The Hi-Loo Angled model is specifically recommended for individuals with restricted knee or hip mobility. Its angled design aims to provide enhanced comfort and accessibility, making it easier for individuals with limited mobility to use the toilet with greater ease and convenience.
- The ergonomic seat of the Hi-Loo Angled model is thoughtfully designed to provide easy access

and comfort. It features recesses at both the back and front of the seat, allowing for convenient and comfortable seating. These recesses enhance accessibility and promote a more pleasant experience for users, particularly those with restricted mobility or specific needs.

## SPECIFICATIONS

Item Name	Item number	Arm support height from top of seat (mm)	CC between bolt holes in toilet bowl (mm)	Front height (mm)	Material	Maximum user weight (kg)	Opening: depth x width (mm)	Rear seat height (mm)	Weight (kg)	Width (mm)	Width between arm supports (mm)
Hi-Loo with arm supports, 6 cm	80301316-2	230	128-202	60	Seat unit: polypropylene. Arm supports: Polypropylene, powder coated steel tubing. Mountings: aluminium.	150	270 x 195	60	4,9	600	490
Hi-Loo with arm supports, 10 cm	80301317-2	230	128-202	100	Seat unit: polypropylene. Arm supports: Polypropylene, powder coated steel tubing. Mountings: aluminium.	150	270 x 195	100	5,8	600	490
Hi-Loo angled with arm supports, 6 cm	80301318-2	230	128-202	60	Seat unit: polypropylene. Arm supports: Polypropylene, powder coated steel tubing. Mountings: aluminium.	150	270 x 195	100	5,5	600	490

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)