




ETAC MY-LOO SOFT SEAT FOR TOILET SEAT RAISER



SKU: 80301604

PRODUCT DESCRIPTION

Purpose & Benefits

The Etac My-Loo Soft Seat is a comfort-enhancing accessory designed for use with the Etac My-Loo fixed toilet seat raiser. Made from soft, cushioned foam with a hygienic, sealed surface, it offers pressure relief and improved comfort for users who require extended time seated on the toilet, ideal for those with sensitive skin, reduced mobility, or in post-operative care.

The seat fits securely onto the My-Loo raiser and is easy to clean, making it a practical addition for aged care, rehabilitation, and home care environments where comfort and hygiene are equally important.

Features

- Soft foam seat for added user comfort and pressure relief
- Designed to fit securely on Etac My-Loo fixed seat raisers
- Hygienic, sealed surface that is easy to wipe clean
- Provides additional comfort during extended toileting
- Simple to fit and remove

- Ideal for elderly users or post-surgical recovery
- Suitable for use in homes, aged care, and clinical settings

Specifications

- **Colour:** Grey
- **Material:** Polyurethane
- **Safe Working Load:** 150 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au