



## **ETAC MY-LOO TOILET SEAT RAISER WITH ARM SUPPORTS**

**\$520.00 - \$630.00** Price range: \$520.00 through \$630.00



**SKU:** 80301511-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Etac My-Loo Toilet Seat Raiser with Arm Supports is a safe, hygienic, and ergonomically designed solution for individuals who require additional height and support when using the toilet. With integrated, softly contoured armrests and three height options, it helps reduce strain on the knees and hips during transfers—making it ideal for elderly users, post-surgical recovery, or individuals with reduced mobility.

The seat is installed using a secure top-mount (blind fixing) system, allowing for tool-free setup on most standard toilet pans. Its smooth surface is easy to clean, and the Scandinavian design ensures it blends discreetly into any bathroom.

### **Features**

- Raises toilet height to assist with sitting and standing
- Integrated ergonomic arm supports for enhanced safety
- Soft contours and wide seating area for added comfort
- Tool-free, top-mounted blind fixing system
- Stable, corrosion-resistant mounting brackets

- Smooth, easy-to-clean surface
- Three height options to suit user needs
- Ideal for aged care, home modifications, and recovery support

## **Specifications**

- Available Heights: 20 mm, 60 mm, 100 mm
- Width Between Arm Supports: 48 cm
- Armrest Height from Seat: 18 cm
- Weight Capacity: 150 kg
- Mount Type: Top-mounted blind fitting
- Materials: Polypropylene seat, powder-coated steel arm supports
- Colour: White