



ETAC MY-LOO TOILET SEAT RAISER WITH LID AND BRACKETS - 60MM

\$190.00



SKU: 80301522

PRODUCT DESCRIPTION

Purpose & Benefits

The **Etac My-Loo Toilet Seat Raiser** with lid and brackets elevates standard toilet seats for enhanced accessibility and comfort. Its ergonomic design, featuring a generous opening and smooth contours, supports safe, independent toileting and seamlessly blends into modern bathrooms.

Features

- Quick, tool-free installation: click-in brackets adjust flexibly to fit most toilet bowls
- Smooth, easy-to-clean surfaces—seat lifts off effortlessly for thorough sanitisation
- Generous front opening and rear recess enhance hygiene and user comfort

Specifications

- **Seat Height:** 6 cm
- **Outer Dimensions (W × D):** 360 mm × 400 mm
- **Opening (D × W):** 360 mm × 230 mm
- **Bracket Range:** 290–380 mm

- **Weight (with lid):** 1.2 kg
- **Maximum User Weight:** 190 kg
- **Materials:** Polypropylene and TPE