



## ETAC RELAX SHOWER SEAT



**SKU:** 8170-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Relax shower seat combines style and functionality. It features a foldable design that saves space by neatly folding it flat against the wall.

Designed with Scandinavian minimalism in mind, the Etac Relax shower seat offers maximum comfort and adds a touch of elegance to your bathroom.

You can customize the seat to meet your specific needs by choosing options such as arm supports, supporting legs, or heat-reflecting back support. The seat itself is designed with a velvety surface and a mild contour, ensuring stability and comfort. Its unique front design allows for alternative sitting positions.

The arm supports have a non-slip coating and ergonomic design, providing a firm grip, while the optional supporting legs offer additional stability if required.

The Etac Relax shower seat is available in white and volcano grey, allowing you to select the colour that suits your bathroom decor.

### FEATURES

- The arm supports of the Etac Relax shower seat are designed with an easy-grip feature. Their ergonomic shape is carefully contoured to fit the hand, ensuring a comfortable and secure grip. This design element enhances safety and stability, allowing users to feel confident and supported while using the shower seat.
- The rounded design of the seat on the Etac Relax shower seat offers multiple sitting positions

and easy access from various sides.

- The Etac Relax shower seat combines stability and a discreet design. It provides a secure and spacious seating surface during use, and when not in use, it folds flat against the wall, maintaining a subtle and space-saving presence in the bathroom.
- The velvety surface of the Etac Relax shower seat provides a smooth and comfortable sitting experience. It is designed to be non-slip, even when wet, ensuring enhanced safety during use.

## SPECIFICATIONS

Item Name	Item number	Colour	Depth (mm)	Folded depth (mm)	Maximum user weight (kg)	Seat height with supporting legs (mm)	Width (mm)	Width outer, with arm supports (mm)
Seat	81703010	White	450	80	125	-	400	-
Seat with arm supports	81703020	White	450	80	125	-	400	530
Seat with supporting legs	81703030	White	450	80	150	450-550	400	-
Seat with arm supports and back support	81703060	White	450	80	125	-	400	530
Seat with arm supports and supporting legs	81703040	White	450	80	150	450-550	400	530
Seat with arm supports, back support and supporting legs	81703050	White	450	80	150	450-550	400	530
Seat	81708000	Volcano grey	450	80	125	-	400	-
Seat with arm supports	81708010	Volcano grey	450	80	125	-	400	530
Seat with supporting legs	81708020	Volcano grey	450	80	150	450-550	400	-
Seat with arm supports and back support	81708050	Volcano grey	450	80	125	-	400	530
Seat with arm supports and supporting legs	81708030	Volcano grey	450	80	150	450-550	400	530
Seat with arm supports, back support and supporting legs	81708040	Volcano grey	450	80	150	450-550	400	530

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)