

ETAC RELIEVE ANGLED CARVING KNIFE

\$43.00



SKU: 80501101



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Relieve angled carving knife is specifically designed to assist individuals with reduced mobility in cutting tasks while standing. The angled design of the knife maximizes strength and efficiency during cutting in an upright position.

The handle of the Relieve carving knife is ergonomically shaped to keep the wrist in a straight position and ensure a stable grip. The top of the handle is smoothly curved to allow for a comfortable thumb position.

To enhance cutting performance, the blade of the knife is sharp and of sufficient length to handle various types of food, including bread. To minimize the risk of contact between the hand and the food, the blade is positioned lower than the handle.

The Etac Relieve angled carving knife provides a practical solution for individuals with reduced mobility, enabling them to perform cutting tasks comfortably and effectively.

FEATURES

- Angled design for optimal strength and efficiency while cutting standing up
- Ergonomically designed handle to keep the wrist straight and ensure a stable grip
- The smoothly curved top of the handle for a comfortable thumb position
- Sharp blade for easy and effective cutting
- Blade is positioned lower than the handle to minimize the risk of hand-food contact
- Sufficient blade length for cutting various types of food, including bread
- Specifically designed for individuals with reduced mobility

- Facilitates cutting tasks while standing, providing increased independence and convenience

SPECIFICATIONS

Item Name	Item number	Blade material	Care	Design	Dimensions	Handle material	Weight (g)
Etac Relieve angled carving knife	80501101	Stainless steel	Dishwasher-safe	Ergonomidesign	Blade: 147 mm. Handle: 109 mm.	Polypropylene	75