

ETAC RELIEVE ANGLED TABLE KNIFE



SKU: 8040210-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Relieve angled table knife assists individuals with limited hand strength in cutting tasks. These knives feature angled handles that make cutting easier for those with reduced hand strength.

You can choose from three variants of the Relieve table knives: one for normal-sized hands, one for small or children's hands, and a foldable model for added convenience. The angled handle design promotes a straight wrist and ensures a stable grip. The top of the handle is smoothly curved to provide a comfortable thumb position.

The blade of the knife is very sharp, facilitating easy cutting. To enhance safety and minimize the risk of hand-food contact, the blade is positioned lower than the handle.

The Etac Relieve angled table knife is a practical solution that enables individuals with limited hand strength to cut food more easily and independently.

FEATURES

- Etac Relieve angled table knife is designed for individuals with limited hand strength
- Angled handles for easier cutting with reduced hand strength
- Available in three variants: normal-sized hands, small/children's hands, and a foldable model
- Angled handle design promotes a straight wrist and stable grip
- The smoothly curved top of the handle for a comfortable thumb position
- Very sharp blade for effortless cutting
- The blade is positioned lower than the handle to minimize hand-food contact
- A practical solution for easier and independent food cutting for individuals with limited hand

strength

SPECIFICATIONS

Item Name	Item number	Blade material	Care	Design	Dimensions	Handle material	Weight (g)
Standard	80402101	Stainless steel	Dishwasher-safe	Ergonomic design	Blade: 117 mm. Handle: 109 mm.	Polyamide plastic	72
Small	80402102	Stainless steel	Dishwasher-safe	Ergonomic design	Blade: 88 mm. Handle: 92 mm.	Polyamide plastic	37

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au