


## ETAC SMART SHOWER STOOL



**SKU:** 819510-H



## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Etac Smart shower stool is a versatile seating solution designed for comfort and stability in the shower. Its rectangular shape offers different seating options, with one side providing a wider seating area and the other side offering a narrower option. The soft corners of the stool make it easy to straddle, ensuring extra stability.

The stool is available in variants including a Swivel pad option and a low version suitable for use as a bath seat or in care situations. With its flexibility, comfort, and practicality, the Etac Smart shower stool is an excellent choice to enhance the showering experience and accommodate various needs.

For a more comfortable sitting experience, we recommend exploring the [Etac Smart Stool Pad](#).

### FEATURES

- Comfortable and safe seating with generous areas and rounded shapes
- Non-slip surface for enhanced safety and comfortable sitting
- Special perforations for quick drainage of excess water
- Height adjustable with easy click-and-press mechanism
- The oval shape prevents telescopic legs from turning during height adjustment
- No loose parts for assembly, tool-free setup
- Self-locking button for secure assembly
- Angled legs for increased stability
- Responsive non-slip ferrules adapt to the floor for stability on wet and uneven surfaces

## SPECIFICATIONS

Item Name	Item number	Colour	Material	Maximum user weight (kg)	Seat area (mm)	Seat height (mm)	Weight (kg)	Width (mm)
Stool	81951010	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	430 x 380	420 - 570	3,2	455 - 465 x 410 - 430
Stool	81951060	White	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	430 x 380	420 - 570	3,2	455 - 465 x 410 - 430
Stool including swivel pad	81951050	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	430 x 380	420 - 570	3,2	455 - 465 x 410 - 430
Stool, low	81951030	Grey	Seat: Polypropylene. Legs: Aluminium, polyamide. Ferrules: TPE	200	430 x 380	345 - 420	3,0	455 - 460 x 410 - 415

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)