

ETAC STAR STABIL-AIR AIR CELL CUSHION



SKU: ESSA-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Star StabilAir cushion offers a unique combination of stability and comfort, making it an excellent choice for those seeking the benefits of an air cushion with the familiar feel of foam.

With its single-chamber vertical air cell design, the StabilAir cushion incorporates foam baffles that contribute to added stability. The foam cylinders within the air cells are strategically positioned, with higher elevation on the outer areas and lower elevation under the buttocks, coccyx, and legs. This configuration slows down the airflow within the cushion, enhancing stability and reducing excessive movement. The air pocket located above the foam also aids in effective pressure redistribution, ensuring optimal comfort and support, as seen in all Star Cushions.

For users who prefer the feel of a foam cushion but require the protection and benefits of a vertical air cell cushion, the StabilAir is the ideal solution. It offers a seamless transition by providing a similar sensation to a foam cushion while delivering the advantages of an air cushion.

The Star StabilAir cushion combines the best of both worlds, offering stability, pressure redistribution, and the familiar feel of a foam cushion, making it an excellent choice for individuals seeking a reliable and comfortable seating solution.

FEATURES

- Vertical Air Cell design for effective pressure redistribution
- Interconnected air cells with high envelopment and immersive properties
- Reduced surface tension for improved comfort
- Excellent pressure redistribution properties

- Continual readjustment to the user's shape as they move in their seat
- Adjustability for optimal immersion with adjustable air levels
- Foam baffles in each cell to slow air movement and increase stability
- Higher foam cylinders on the side for enhanced stability

SPECIFICATIONS

Item number	Cell height (mm)	Maximum user weight	Suits wheelchair seat (mm)	
26 x 26 cm	ESSA30909-1	75	Unlimited	220-250
26 x 30 cm	ESSA30911-1	75	Unlimited	220-250
26 x 34 cm	ESSA30913-1	75	Unlimited	220-250
30 x 26 cm	ESSA31109-1	75	Unlimited	250-300
30 x 30 cm	ESSA31111-1	75	Unlimited	250-300
30 x 34 cm	ESSA31113-1	75	Unlimited	250-300
30 x 39 cm	ESSA31115-1	75	Unlimited	250-300
34 x 26 cm	ESSA31309-1	75	Unlimited	300-350
34 x 30 cm	ESSA31311-1	75	Unlimited	300-350
34 x 34 cm	ESSA31313-1	75	Unlimited	300-350
34 x 39 cm	ESSA31315-1	75	Unlimited	300-350
34 x 44 cm	ESSA31317-1	75	Unlimited	300-350
39 x 30 cm	ESSA31511-1	75	Unlimited	350-400
39 x 34 cm	ESSA31513-1	75	Unlimited	350-400
39 x 39 cm	ESSA31515-1	75	Unlimited	350-400
39 x 44 cm	ESSA31517-1	75	Unlimited	350-400
39 x 49 cm	ESSA31519-1	75	Unlimited	350-400
39 x 53 cm	ESSA31521-1	75	Unlimited	350-400
44 x 34 cm	ESSA31713-1	75	Unlimited	400-440
44 x 39 cm	ESSA31715-1	75	Unlimited	400-440
44 x 44 cm	ESSA31717-1	75	Unlimited	400-440
44 x 49 cm	ESSA31719-1	75	Unlimited	400-440
44 x 53 cm	ESSA31721-1	75	Unlimited	400-440
49 x 39 cm	ESSA31915-1	75	Unlimited	440-490
49 x 44 cm	ESSA31917-1	75	Unlimited	440-490
49 x 49 cm	ESSA31919-1	75	Unlimited	440-490
49 x 53 cm	ESSA31921-1	75	Unlimited	440-490
53 x 39 cm	ESSA32115-1	75	Unlimited	490-540

53 x 44 cm	ESSA32117-1	75	Unlimited	490-540
53 x 49 cm	ESSA32119-1	75	Unlimited	490-540
53 x 53 cm	ESSA32121-1	75	Unlimited	490-540
26 x 26 cm	ESSA40909-1	100	Unlimited	220-250
26 x 30 cm	ESSA40911-1	100	Unlimited	220-250
26 x 34 cm	ESSA40913-1	100	Unlimited	220-250
30 x 26 cm	ESSA41109-1	100	Unlimited	250-300
30 x 30 cm	ESSA41111-1	100	Unlimited	250-300
30 x 34 cm	ESSA41113-1	100	Unlimited	250-300
30 x 39 cm	ESSA41115-1	100	Unlimited	250-300
34 x 26 cm	ESSA41309-1	100	Unlimited	350-400
34 x 30 cm	ESSA41311-1	100	Unlimited	300-350
34 x 34 cm	ESSA41313-1	100	Unlimited	300-350
34 x 39 cm	ESSA41315-1	100	Unlimited	300-350
34 x 44 cm	ESSA41317-1	100	Unlimited	300-350
39 x 30 cm	ESSA41511-1	100	Unlimited	350-400
39 x 34 cm	ESSA41513-1	100	Unlimited	350-400
39 x 39 cm	ESSA41515-1	100	Unlimited	350-400
39 x 44 cm	ESSA41517-1	100	Unlimited	350-400
39 x 49 cm	ESSA41519-1	100	Unlimited	350x400
39 x 53 cm	ESSA41521-1	100	Unlimited	350-400
44 x 34 cm	ESSA41713-1	100	Unlimited	400-440
44 x 39 cm	ESSA41715-1	100	Unlimited	400-440
44 x 44 cm	ESSA41717-1	100	Unlimited	400-440
44 x 49 cm	ESSA41719-1	100	Unlimited	400-440
44 x 53 cm	ESSA41721-1	100	Unlimited	400-400
49 x 39 cm	ESSA41915-1	100	Unlimited	440-490
49 x 44 cm	ESSA41917-1	100	Unlimited	440-490
49 x 49 cm	ESSA41919-1	100	Unlimited	440-490
49 x 53 cm	ESSA41921-1	100	Unlimited	440-490
53 x 39 cm	ESSA42115-1	100	Unlimited	490-540
53 x 44 cm	ESSA42117-1	100	Unlimited	490-540
53 x 49 cm	ESSA42119-1	100	Unlimited	490-540
53 x 53 cm	ESSA42121-1	100	Unlimited	490-540

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au