

ETAC STAR STANDARD AIR CELL CUSHION



SKU: ESSC-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Standard Air cushion is a traditional single-chamber vertical air cell cushion designed to provide pressure redistribution capabilities.

The Standard Air cushion offers adjustability to ensure optimal immersion for the user, resulting in excellent pressure redistribution. Its interconnected air cell structure allows for unrestricted airflow within the cushion, allowing it to continually adapt to the user's body shape as they move in their seat. This dynamic feature promotes enhanced comfort and support throughout prolonged sitting periods.

FEATURES

- The interconnected design of the vertical air cells allows for effective weight distribution and pressure relief, resulting in improved comfort and a decreased risk of pressure injuries.
- The air level of the cushion can be adjusted and reset repeatedly to achieve optimal immersion for the individual user. This flexibility allows for personalized customization, ensuring comfort and support based on individual needs.
- The interconnected air cells in the cushion continually readjust to the user's shape as they move, ensuring excellent pressure redistribution properties.
- The Standard Air Cushion features a single chamber design, which means it requires only one valve adjustment. This design simplifies the handling and maintenance of the cushion, making it ideal for easy use and convenience.

SPECIFICATIONS

Item Variant	Item number	Cell height (mm)	Maximum user weight	Suits wheelchair seat (mm)
26 x 26 cm	ESSCS0909-1	50	Unlimited	220-250
26 x 30 cm	ESSCS0911-1	50	Unlimited	220-250
26 x 34 cm	ESSCS0913-1	50	Unlimited	220-250
30 x 26 cm	ESSCS1109-1	50	Unlimited	250-300
30 x 30 cm	ESSCS1111-1	50	Unlimited	250-300
30 x 34 cm	ESSCS1113-1	50	Unlimited	250-300
30 x 39 cm	ESSCS1115-1	50	Unlimited	250-300
34 x 26 cm	ESSCS1309-1	50	Unlimited	300-350
34 x 30 cm	ESSCS1311-1	50	Unlimited	300-350
34 x 34 cm	ESSCS1313-1	50	Unlimited	300-350
34 x 39 cm	ESSCS1315-1	50	Unlimited	300-350
34 x 44 cm	ESSCS1317-1	50	Unlimited	300-350
39 x 30 cm	ESSCS1511-1	50	Unlimited	350-400
39 x 34 cm	ESSCS1513-1	50	Unlimited	350-400
39 x 39 cm	ESSCS1515-1	50	Unlimited	350-400
39 x 44 cm	ESSCS1517-1	50	Unlimited	350-400
39 x 49 cm	ESSCS1519-1	50	Unlimited	350-400
39 x 53 cm	ESSCS1521-1	50	Unlimited	350-400
44 x 34 cm	ESSCS1713-1	50	Unlimited	400-440
44 x 39 cm	ESSCS1715-1	50	Unlimited	400-440
44 x 44 cm	ESSCS1717-1	50	Unlimited	400-440
44 x 49 cm	ESSCS1719-1	50	Unlimited	400-440
44 x 53 cm	ESSCS1721-1	50	Unlimited	400-440
49 x 39 cm	ESSCS1915-1	50	Unlimited	440-490
49 x 49 cm	ESSCS1919-1	50	Unlimited	440-490
49 x 53 cm	ESSCS1921-1	50	Unlimited	440-490
49 x 44 cm	ESSCS1917-1	50	Unlimited	440-490
53 x 44 cm	ESSCS2117-1	50	Unlimited	490-540
53 x 39 cm	ESSCS2115-1	50	Unlimited	490-540
53 x 49 cm	ESSCS2119-1	50	Unlimited	490-540
53 x 53 cm	ESSCS2121-1	50	Unlimited	490-540
26 x 26 cm	ESSC0909-1	100	Unlimited	220-250
26 x 30 cm	ESSC0911-1	100	Unlimited	220-250

26 x 34 cm	ESSC0913-1	100	Unlimited	220-250
30 x 26 cm	ESSC1109-1	100	Unlimited	250-300
30 x 30 cm	ESSC1111-1	100	Unlimited	250-300
30 x 34 cm	ESSC1113-1	100	Unlimited	250-300
30 x 39 cm	ESSC1115-1	100	Unlimited	250-300
34 x 26 cm	ESSC1309-1	100	Unlimited	300-350
34 x 30 cm	ESSC1311-1	100	Unlimited	300-350
34 x 34 cm	ESSC1313-1	100	Unlimited	300-350
34 x 39 cm	ESSC1315-1	100	Unlimited	300-350
34 x 44 cm	ESSC1317-1	100	Unlimited	300-350
39 x 30 cm	ESSC1511-1	100	Unlimited	350-400
39 x 34 cm	ESSC1513-1	100	Unlimited	350-400
39 x 39 cm	ESSC1515-1	100	Unlimited	350-400
39 x 44 cm	ESSC1517-1	100	Unlimited	350-400
39 x 49 cm	ESSC1519-1	100	Unlimited	350-400
39 x 53 cm	ESSC1521-1	100	Unlimited	350-400
44 x 34 cm	ESSC1713-1	100	Unlimited	400-440
44 x 39 cm	ESSC1715-1	100	Unlimited	400-440
44 x 44 cm	ESSC1717-1	100	Unlimited	400-440
44 x 49 cm	ESSC1719-1	100	Unlimited	400-440
44 x 53 cm	ESSC1721-1	100	Unlimited	400-440
49 x 39 cm	ESSC1915-1	100	Unlimited	440-490
49 x 44 cm	ESSC1917-1	100	Unlimited	440-490
49 x 49 cm	ESSC1919-1	100	Unlimited	440-490
49 x 53 cm	ESSC1921-1	100	Unlimited	400-440
53 x 39 cm	ESSC2115-1	100	Unlimited	490-540
53 x 44 cm	ESSC2117-1	100	Unlimited	490-540
53 x 49 cm	ESSC2119-1	100	Unlimited	490-540
53 x 53 cm	ESSC2121-1	100	Unlimited	490-540

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au