



ETAC SUPPORTER TOILET ARMS

\$558.00



SKU: 80303118-2

PRODUCT DESCRIPTION

Purpose & Benefits

The Etac Supporter Toilet Arms offer discreet, stable assistance for individuals who need extra support when using the toilet. Designed with comfort, safety, and minimal intrusion in mind, these armrests provide a practical solution for anyone with reduced strength, balance issues, or recovering from surgery.

Featuring flip-up arms and adjustable support legs, the Etac Supporter promotes easier, safer transfers in and out of the seated position. With clean Scandinavian styling and tool-free installation on most standard toilet bowls, it's an ideal choice for aged care, home modifications, and accessible bathroom upgrades.

Features

- Provides secure arm support for toilet use
- Flip-up armrests for easy side transfers and cleaning access
- Height-adjustable support legs for extra stability
- Ergonomic, softly contoured armrest grips
- Durable, rust-resistant materials for long-term use

- Compact design blends into modern bathrooms
- Easy to install—fits most standard toilets
- Ideal for aged care and home safety solutions

Specifications

- Width Between Arm Supports: 48 cm
- Armrest Height from Seat: 26 cm
- Safe Working Load: 150 kg
- Adjustment: Height-adjustable legs
- Armrest Function: Flip-up for access
- Material:
 - Frame and arm supports: powder coated tubular steel.
 - Arm supports cover: polypropylene.
 - Lid and seat: Polypropylene.
 - Mountings: Aluminium
- Installation: Mounts directly to toilet bowl