



## ETAC TOILET SEAT SUPPORTER WITH ARM SUPPORTS



**SKU:** 80303112-2

## PRODUCT DESCRIPTION

### Purpose & Benefits

The Etac Supporter is a discreet and supportive toilet armrest system designed to improve safety and independence in the bathroom. Ideal for individuals with reduced mobility, balance issues, or post-operative needs, the Supporter provides sturdy, ergonomic support when sitting down or standing up from the toilet.

Its flip-up armrests offer convenience for transfers and space-saving flexibility, while the height-adjustable legs provide stable, customisable support. With clean Scandinavian design and easy installation on most standard toilet bowls, the Supporter is a practical and non-intrusive solution for home or aged care use.

### Features

- Provides sturdy support for toilet transfers
- Flip-up armrests for easy side access and compact storage
- Height-adjustable legs for added stability and custom fit
- Non-slip armrest surfaces for secure grip
- Clean, modern design fits most bathrooms discreetly

- Easy to install and remove without major modifications
- Ideal for elderly users or those with mobility challenges

## Specifications

- Armrest Height from Seat: 260 mm
- Width Between Armrests: 480 mm
- Weight Capacity: 130 kg
- Armrest Flip-Up Function: Yes
- Height Adjustment (legs): Yes
- Material: Powder-coated steel with polyurethane armrests
- Compatibility: Fits most standard toilets

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)