



ETAC TOILET SEAT SUPPORTER WITH ARM SUPPORTS

\$480.00



SKU: 80303112-2

PRODUCT DESCRIPTION

Purpose & Benefits

The Etac Supporter is a discreet and supportive toilet armrest system designed to improve safety and independence in the bathroom. Ideal for individuals with reduced mobility, balance issues, or post-operative needs, the Supporter provides sturdy, ergonomic support when sitting down or standing up from the toilet.

Its flip-up armrests offer convenience for transfers and space-saving flexibility, while the height-adjustable legs provide stable, customisable support. With clean Scandinavian design and easy installation on most standard toilet bowls, the Supporter is a practical and non-intrusive solution for home or aged care use.

Features

- Provides sturdy support for toilet transfers
- Flip-up armrests for easy side access and compact storage
- Height-adjustable legs for added stability and custom fit
- Non-slip armrest surfaces for secure grip
- Clean, modern design fits most bathrooms discreetly

- Easy to install and remove without major modifications
- Ideal for elderly users or those with mobility challenges

Specifications

- Armrest Height from Seat: 260 mm
- Width Between Armrests: 480 mm
- Weight Capacity: 130 kg
- Armrest Flip-Up Function: Yes
- Height Adjustment (legs): Yes
- Material: Powder-coated steel with polyurethane armrests
- Compatibility: Fits most standard toilets