

ETAC SWIFT MOBIL TILT-2 SHOWER COMMODE



SKU: 802294-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Engineered keeping in mind the needs of both caregivers and users, this product addresses the vast array of challenges caregivers encounter daily while ensuring user comfort and safety.

Exclusive Tilt-in-Space Feature - The Swift Mobil Tilt-2 possesses a distinctive tilt mechanism that enhances the user's security and provides an ergonomic position for the caregiver. Its -5° to $+35^{\circ}$ tilt functionality can be controlled from either side of the chair, allowing the caregiver to comfort the user during tilting.

Accessibility and Comfort - Outstanding access around the seat area simplifies personal care routines for the caregiver while maintaining the user's dignity. Extra-soft seat cushions are available for sensitive users. The soft, fast-drying backrest and curved footrests of the Swift Mobil-2 contribute to the user's comfort.

Suitable for a Variety of Users - The shower commode chair is versatile and can be adjusted to accommodate a wide range of user needs and sizes. It supports a generous maximum user weight of 160 kg as a standard feature and offers adjustable seat-to-floor height (50-65 cm). The armrests can be adjusted in terms of both height and width, and the height-adjustable, swing-away footrests are also available in a shorter version.

Etac Swift Mobil Tilt-2 XL for Larger Users - A variant with a broader back frame and leg rest widening adapters is available, offering an extra 6 cm of space between armrests and an additional 12 cm of legroom.

FEATURES

- Seat with rear opening - Ensuring convenient accessibility for the caregiver and comfort for the user. The one-of-a-kind surface, specifically created for Etac, strikes the ideal balance between grip and cleanliness, offering both non-slip attributes and hygienic conditions!
- Adjustable armrests as standard - Effortlessly customizable to meet diverse user requirements. It offers both height and width adjustments, with a range of +/- 6 cm in each direction.
- Unique and safe tilting function - The chair can be controlled from either side, facilitating an ergonomic working posture for the caregiver and providing comfort to the user during the tilting process.
- Multi-adjustable headrest - The headrest is meticulously designed to ensure optimal comfort and versatility. It can be adjusted to accommodate a wide range of user needs, allowing for a personalized and comfortable experience.

SPECIFICATIONS

| Item Name | Item number | Length (mm) | Material | Maximum user weight (kg) | Opening: depth x width (mm) | Seat depth (mm) | Seat height (mm) | Seat width (mm) | Space below seat (mm) | Tilt in space tilt angle (°) | Weight (kg) | Width (mm) | Width between arm supports (mm) |
|--|-------------|-------------|---|--------------------------|-----------------------------|-----------------|------------------|-----------------|-----------------------|--|-------------|-------------|---------------------------------|
| Standard Swift Mobil Tilt-2, with pan holder | 80229412 | 1020 | Chassis: powder-painted steel. Seat/Arm supports: polypropylene. Back support: polyester/velcro. Foot support: polypropylene. Castors: polyamide/TPE. | 160 | 300x210 | 465 | 500-650 | 500 | 427-575 | 0° t/m +35° /-5° t/m +30° | 20,4 | 570/580/640 | 420/480/540 |
| Swift Mobil Tilt-2 XL-back, with pan holder | 80229481 | 1020 | Chassis: powder-painted steel. Seat/Arm supports: polypropylene. Back support: polyester/velcro. Foot support: polypropylene. Castors: polyamide/TPE. | 160 | 300 x 210 | 465 | 500 - 650 | 500 | 427 - 575 | 0° t/m +35° /-5° t/m +30° | 21 | 610/640/700 | 480/540/600 |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au