


ETAC SWIFT FREE STANDING TOILET SEAT RAISER



SKU: 81702020



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Swift Freestanding is the optimal selection for a raised toilet seat that can be effortlessly set up. It provides a customizable seat height that caters to your individual requirements and robust armrests that simplify the motions of standing and sitting. When augmented with a pan, it can also function in the bedroom as an independent toilet chair. Due to its design, which is resistant to damp environments, the Swift can double as a shower stool.

FEATURES

- **Compatibility with diverse toilets:** Swift can be seamlessly adjusted to the required height and accommodates both wall and floor-mounted toilets. The sturdy armrests render significant support for the user during the process of standing up or sitting down.
- **Impressive stability:** The chair's broad footprint contributes to its exceptional stability. The adaptable, non-slip ferrules on its legs provide an excellent grip on wet or irregular surfaces.
- **Effortless setup:** The assembly process of the chair is straightforward, requiring no tools. The self-locking mechanism snaps into place securely, and the oval-shaped design of the legs avoids any rotation during adjustments.
- **Easy transportation:** Swift is packaged in a flat, compact carry case for hassle-free transit. Given its ease of assembly and disassembly, the case can be retained and utilized for additional transports if the chair is used in various locations.

SPECIFICATIONS

Seat height	42–57 cm
--------------------	----------

Seat width	54 cm
Width between arm supports	45 cm
Width	56 cm
Opening: depth x width	30 cm x 21,5 cm
Weight	3.9 kg
Maximum user weight	130 kg
Material	Seat and arm supports: polypropylene Legs: Aluminum and polyamide Ferrules: TPE

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au