



ETAC TURNER PRO

\$1,125.00



SKU: 16090106

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Turner PRO is the perfect tool for hospitals, aged care facilities or home care use, for patients who are able to bear their own weight when standing. The intuitive design makes it easy to use, the only adjustments required being handle height and height of leg supports.

Reduces risk of patient falls and ensures carers can maintain an ergonomic working position. Patient actively moves to a standing position, promoting independence and continued mobility. Allows carers to work ergonomically and maintain eye contact, while patient controls movement for enhanced security.

FEATURES

- Padded, Indented Knee Pads for Added Safety During Transfer
- Wheels Offer Easier Transport Across Different Floor Types
- Quickly and Easily Rises a Patient to a Standing Position
- Great for Transferring from Bed to Chair
- Swivel Bottom Assures Quick Transfer of the Patient
- Dual Handle Bar Offers Room for Both Parties
- Lightweight for Easier Carrying if Needed
- Bright-Colored Handle for Easy Viewing