



EXOFORM CARPAL TUNNEL WRIST

\$39.00



SKU: OSS517073-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Exform Carpal Tunnel Wrist is a lightweight and low profile wrist brace that provides immobilisation to the joint.

FEATURES

- Lightweight low-profile design
- More room at base of the thumb, palmar crease and 5th MP for full MP flexion
- Designed to provide smooth controlled resistance to extreme range of motion
- Pre-shaped into a neutral position for comfortable limited motion/flexibility
- Provides radial and ulnar styloid relief
- Soft touch web space is cushioned and 40% thinner than other brands
- Interior is 3-dimensionally moulded to enhance comfort

SPECIFICATIONS

| Wrist Circumference* | Palm Circumference** | |
|---------------------------------|---------------------------------|-----------|
| Small | 13 - 16.5cm | 17 - 20cm |
| Medium | 16.5 - 20cm | 20 - 23cm |
| Large | 20 - 25cm | 23 - 27cm |