



MERITS EZI-GO DLX



SKU: P321A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Merits Ezi-Go Deluxe Powerchair dismantles in seconds for easy transport and storage. The tyres are all puncture proof and the chair is suitable for indoor and shopping centre use. The armrests are height and width adjustable and the brakes are activated automatically when the joy stick is released. The Ezi-Go Deluxe has a maximum speed of up to 6.4km/hr and a range of up to 15km.

FEATURES

- Battery and charger included
- very manoeuvrable, ideal for indoor use
- Super light and easy disassembly for carrying, storage and transportation
- Compact size to easily fit into the boot of most cars
- High strength steel frame
- Comfortable seat with adjustable height & width armrest
- Quick release trigger separates the base sections
- Colours available: Red or Blue
- DLX (low-back) and XL (high-back) sizes available