


FINCH METZELER PRESSURE CARE MATTRESS



SKU: FM-2AD-076200-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Metzeler Pressure Care Mattress

The Metzeler Pressure Care Mattress is meticulously designed to provide superior comfort, support, and effective pressure relief for individuals requiring enhanced care during sleep. Incorporating innovative Rubex™ Duo Cell technology, this mattress adapts to the body's contours, ensuring optimal pressure distribution and promoting spinal alignment. Ideal for those at risk of pressure ulcers or seeking improved sleep quality, the Metzeler mattress offers a restorative sleep experience.

Features

- **Rubex™ Duo Cell Technology:** Patented dual-layer cell structure combines large and small pores to reduce pressure and provide soft support, enhancing air permeability and moisture transport for consistent bed temperature.
- **Multizone Support System:** Unique tubular core structure zoned to support the spine, offering optimal pressure redistribution and reducing stress on pressure points.
- **Reinforced Mattress Edge:** Durable foam edges assist with safe transfers on and off the bed, minimising the potential for rolling and aiding caregivers during clinical care.
- **Waterproof and Breathable Cover:** Visco-elastic cover aids moisture management, is easy to remove, and washable up to 95°C, ensuring a hygienic sleep surface.
- **Optimal Temperature Management:** Tubular structure and choice of foam overlays provide cross-ventilation horizontally and vertically, facilitating effective thermoregulation.
- **Environmentally Friendly Design:** Manufactured from natural, non-toxic, non-petroleum-based materials, including sunflower oil, making it an eco-friendly alternative to conventional foam mattresses.

Specifications

- **Available Sizes:** Single, King Single, Double, Queen, Queen Split
- **Firmness Options:** Soft, Medium, Firm, Extra Firm
- **Minimum User Weight:** 20 kg
- **Core Height:** 12 cm, with an option for 14 cm
- **Weight Capacity:** Suitable for users up to 180 kg
- **Cover Material:** Air-permeable polyurethane coating, free from latex and PVC
- **Cleaning:** Cover washable up to 95°C; suitable for steam sterilisation up to 105°C

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au