



FISHER & PAYKEL BREVIDA NASAL PILLOW MASK



SKU: BRE1SMA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Fisher & Paykel Brevida is a revolutionary nasal pillows mask that has been designed for confidence and built for comfort. It combines a gentle, soft AirPillow seal which inflated to form a stabilising side support to seal in and around the nose. It features an adjustable headgear to fit all sizes, and a washable diffuser to reduce noise and air draft.

FEATURES

- AirPillow Seal; inflates to form a gentle 'pillow' of air that seals around the nose. Fits a wide variety of noses
- Adjustable headgear; simplistic design that stays in place during sleep
- Washable Diffuser; reduces noise, partner disturbance and ongoing consumable costs
- 2x size cushions included extra small/small and medium/large

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au