



## **FISHER & PAYKEL ESON 2 NASAL MASK**

**\$298.00**



**SKU: ESN2-H**

## **PRODUCT DESCRIPTION**

## **PURPOSE & BENEFITS**

The Eson 2 is the enhanced version of Fisher & Paykel's best performing nasal mask. There are only 3 components to the Eson 2, all designed to work in harmony for greater performance, stability and comfort.

## **FEATURES**

- RollFit Seal; relieves pressure on the bridge of nose while maintaining seal and comfort
- Easy Frame; the light and narrow profile for clearer line of sight
- Intuitive Headgear; straps can be stretched over the head unlike traditional headgear
- Washable Diffuser; reduces noise, partner disturbance and ongoing consumable costs