


FISHER & PAYKEL EVORA NASAL MASK - SMALL

~~\$308.00~~ Original price was: \$308.00. ~~\$223.00~~ Current price is: \$223.00. Sale Available Online Only! While Stocks Last



SKU: EVO1SU



PRODUCT DESCRIPTION

Purpose & Benefits

The **Fisher & Paykel Evora Nasal Mask - Small** is a compact nasal CPAP mask designed for comfort, simplicity, and freedom of movement. Featuring CapFit™ headgear that can be put on like a cap, it reduces complexity and makes setup intuitive. Dynamic Support Technology™ with stability wings and a floating seal helps maintain a secure fit while minimising facial contact, supporting better sleep therapy and comfort.

Features

- CapFit™ headgear designed to be put on like a cap for intuitive fitting
- Dynamic Support Technology™ with stability wings and floating seal for secure comfort
- Compact under-the-nose cushion for minimal contact and a free field of vision
- Exhaust holes engineered to reduce noise and draft

Specifications

- **Size:** Small
- **Available Cushion Sizes:** Small, Medium, Large, Wide

- **Headgear:** One size fits all cushions