




## FISHER & PAYKEL SIMPLUS FULL FACE MASK

**\$308.00**Online Only



**SKU:** 40047-H

## PRODUCT DESCRIPTION

The Fisher & Paykel Simplus Full Face CPAP Mask is the result of intensive research and development, many prototypes and thousands of design hours. Fisher & Paykel consulted patients, equipment providers and clinicians from around the world to develop a mask that simply fits + performs.

**F&P Simplus™, Simply Fits Fast:** Only 3 key components: RollFit Seal, ErgoForm Headgear, Easy Frame. Only 3 adjustments required to effectively fit the F&P Simplus for the first time.

**F&P Simplus™, Simply Cleans:** Quick & Easy: Has Minimal number of parts to clean, Easy reassembly with Easy-Clip Seal, One size frame fits all three seal sizes, (Small, Medium, Large).

All these components are designed to work in harmony. In combination, these components offer the comfort, seal and easy use that Fisher & Paykel Healthcare masks are known for.

## 3 Key Features of the F&P Simplus

**Easy Frame:** Easy to wear, see and sleep. This low-profile frame is stable, durable, small and ensures a clear line of sight. This is a **One-Frame-Fits-All-Seal-Sizes** and has an 'Easy-Clip' Frame attachment for effortless assembly after cleaning. The **Ball & Socket Elbow** rotates for freedom of movement reducing drag on the tube.

**ErgoForm™ Headgear:** Stretch and non-stretch panels provide structure and support for optimal usability and performance. This breathable headgear self-locates high on the rear of the head allowing for maximum head movement (side-to-side, up and down) without mask movement.

**Velcro Tabs** Has **Easy Clip Hooks** that enable fast fitting and removal.

**RollFit™ Seal:** Auto-adjusting to optimize an effective, comfortable seal. The one piece seal 'rolls' back and forth on the bridge of the nose. This automatic adjustment minimizes pressure on the bridge of the nose.