



## FISHER & PAYKEL SLEEPSTYLE FILTERS 2 PACK

\$14.00 Original price was: \$14.00.\(\frac{\pmathbf{\$}}{2.00}\)Current price is: \$12.00.\(\pmathbf{While Stocks LastOnline Only}\)



**SKU:** 900SPS110

## **PRODUCT DESCRIPTION**

## **PURPOSE & USES**

Made from a durable, washable material these filters can last 3-6 months before needing to be replaced.

They are located at the back of the machine and are easily pushed in or pulled out.