



FISHER & PAYKEL UNIVERSAL 19MM ID TUBE (1.8M)

\$34.00 Original price was: \$34.00.\(\frac{\$24.00}{2}\) Current price is: \$24.00.\(\text{While Stocks LastOnline Only}\)



SKU: 900HC221

PRODUCT DESCRIPTION

The 1.8m Long 19mm Diameter CPAP Hose with 22mm Rubber Ends is comparable to the standard, non-heated 6 foot hose that comes with many CPAP and BiLevel machines. The hose measures 19mm interior diameter and a 22mm diameter rubber cuff. The standard cuff size allows the hose to connect to any machine compatible with a standard hose.