


**FISHER & PAYKEL UNIVERSAL 19MM ID TUBE
(1.8M)****\$34.00Online Only****SKU: 900HC221****PRODUCT DESCRIPTION**

The 1.8m Long 19mm Diameter CPAP Hose with 22mm Rubber Ends is comparable to the standard, non-heated 6 foot hose that comes with many CPAP and BiLevel machines. The hose measures 19mm interior diameter and a 22mm diameter rubber cuff. The standard cuff size allows the hose to connect to any machine compatible with a standard hose.