



## FORTE HEALTHCARE ICON MAXX BARIATRIC MATTRESS - KING SINGLE



**SKU:** IM1KS1000

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Forte Healthcare Icon Maxx Bariatric Mattress - King Single** is designed to deliver advanced pressure care, support, and comfort for bariatric patients. Featuring multi-layered foam construction with memory foam immersion, hinged reinforced sides, and postural support layers, it promotes safety, stability, and effective pressure redistribution for long-term care applications.

### Features

- **Airoform Memory Immersion Interface Layer (Medium):** Provides consistent comfort and temperature stability while maximising patient contact for effective pressure redistribution.
- Visco-elastic memory foam structure offers a low-shear surface that alleviates localised pressure points.
- Depth calibrated to balance immersion, functionality, and patient mobility.
- **Hinged Strengthened Sides:** Improve lateral stability and safety, promote central patient positioning, and assist with transfers.
- Unique hinge profile reduces foam fatigue and conforms easily to adjustable bed positions for ease of care.
- **Premium Postural Support Secondary Layer:** Provides graded support and comfort in supine, prone, lateral, and fowler positions.

- Enhances stability to prevent falls, minimise entrapment, and maintain central alignment.
- **Premium Firm Foundation Layer:** Protects against bottoming-out during repositioning and provides consistent support for patient care.

## Specifications

- **Model:** Icon Maxx Bariatric Mattress
- **Size:** King Single
- **Firmness:** Medium (Airoform Memory Immersion Interface Layer)
- **Core Structure:** Multi-layered foam with hinged sides and reinforced base
- **Key Benefits:** Pressure redistribution, postural support, comfort, and safety for bariatric users

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)